



University of Milan
Department of Pharmaceutical Sciences



Chemical Laboratory Simulator of
Qualitative Inorganic Analysis

Alessandro Pedretti

What is LabSim¹




LabSim is a laboratory simulator of qualitative inorganic analysis which was mainly developed with the aim to teach the students in laboratory during the COVID pandemic emergency. Its main features are:

- **Complete 3D environment** (based on A-Frame/WebGL technology).
- **Web-based application:**
 - Compatibility with widest range of devices as: **VR systems** (HTC and Oculus), **PCs** (Linux, MacOS, and Windows), **tablets** (Android, iOS and Windows) and **smartphones** (Android and iOS).
 - Runs inside the web browser (Chromium-based), so **no software installation is required**.



What is LabSim²

- Support multiple languages (actually *English* and *Italian*).
- Speech synthesis (based on different technologies according to the device).
- Speech recognition for disabled persons.
- Based on a database of more than 3500 reactions that can be used in recognizing 60 analytes.



- 70 video virtual and real tutorials watchable through a TV included in the 3D environment.
- Totally free.

What you can do



With *LabSim*, you can practice as in a real laboratory of qualitative inorganic analysis and the main activities are:

- **Solubility test** in water of a substance even in the presence of reactants such as strong / weak acids / bases;
- **pH measurement** of a solution with universal indicator paper;
- **Phase separation** by centrifugation;
- **Heating** by plate or water bath;
- **Tests for the recognition of anions** by wet method;
- **Tests for the recognition of cations** by wet method;
- **Recognition tests by dry method** (non water-soluble substances);
- **Flame test**;
- **Complete analysis of an unknown substance** (both water-soluble and non-water-soluble);
- **Watching video tutorials** on both virtual and real laboratory activities.




How to try it



In the address bar of the browser (Chrome or Edge), type:

[*https://nova.disfarm.unimi.it/labsim*](https://nova.disfarm.unimi.it/labsim)

After the loading screen, it will shown the copyright message:



Click the **Ok** button to start the simulation. Clicking the **VR** button, you switch from widowed mode to full-screen or VR mode. Press **Esc** key to return to the standard display mode.

The mouse




The mouse is the preferred input device for non-VR hardware systems:

Scroller
Zoom in/out the scene

Right button
Keep pressed and move the mouse to translate the scene

Left mouse
Keep pressed and move the mouse to rotate the scene

Click an object to make an action



Other input devices

Some functions are available through the **keyboard**:

- The **cursor key** translate the scene.
- The **space bar** and **m** key show the main menu.
- The **Esc** key closes the menu.




Input devices available only with **tablets** and **smartphones**:


- The **screen touch** (tap) is equivalent to a mouse click.
- The **accelerometer** rotates the scene.



The workbench



The centrifuge



To perform a complete phase separation, **at least 3 minutes** are needed. Watch the video-tutorial for the right use.

How to watch the video tutorials

When you click the TV, the video *playlist* is shown:

Click one of these buttons to play

Click here to turn off the TV

List of the videos

15. LS Riconoscimento PO_4^{3-}

16. LS Riconoscimento CH_3COO^-

17. LS Riconoscimento F^-

18. LS Riconoscimento BO_3^{3-}

19. LS Riconoscimento I gruppo (Ag^+)

20. LS Riconoscimento II gruppo (Hg^{2+})

21. LS Riconoscimento III gruppo (Al^{3+})

22. LS Riconoscimento IV gruppo (Zn^{2+})

23. LS Riconoscimento V gruppo (Ca^{2+})

24. LS Riconoscimento VI gruppo (NH_4^+)

25. LS Riconoscimento VI gruppo (Mg^{2+})

26. LS Riconoscimento VI gruppo (Li^+)

27. LS Riconoscimento VI gruppo (K^+)


28. LS Riconoscimento VI gruppo (Na^+)



Previous page

Next page

How to control the player



Clicking **Stop**, you return to the video playlist.


Direct access to the videos



You can play the video tutorials outside *LabSim* just with this link:

<https://nova.disfarm.unimi.it/labsim/video.htm>

After the loading screen, the following menu is shown:



This playlist works exactly as in *LabSim*.

What you can watch

Tutorials on the virtual lab (LS prefix):

- Base techniques (solubility, pH measurement, centrifuge use).
- Tests for the identifications of the anions.
- Tests for the identifications of the cations.
- Flame test.
- Analysis of water-insoluble molecules.
- Two examples of analysis of unknown molecules.

Tutorials on the real lab:

- Tests for the identifications of the anions.
- Tests for the identifications of the cations.
- Flame test.
- Analysis of some water-insoluble molecules.



How to change the position


Clicking the footprints, you can change the work position.



Main menu



To open it, you must click the book or press *m* key.



Menu analita



In this menu, you can choose the analyte.

Back to the previous menu

Random selection of the analyte

Close the menu

Analyte menu

Code	Random	Random soluble	Random insoluble	$(\text{CH}_3\text{COO})_2\text{Ca}$	$(\text{CH}_3\text{COO})_2\text{Mg}$
AgNO_3	$\text{Al}_2(\text{SO}_4)_3$	AlCl_3	CaCl_2	CH_3COOK	CH_3COONa
H_3BO_3	HgCl_2	K_2CO_3	K_2HPO_4	$\text{K}_2\text{S}_2\text{O}_3$	$\text{K}_2\text{S}_2\text{O}_5$
K_2SO_3	K_2SO_4	KAlSO_4	KBr	KCl	KH_2PO_4
KHCO_3	KI	KIO_3	KNO_3	Li_2CO_3	MgCl_2
MgSO_4	$\text{Na}_2\text{B}_4\text{O}_7$	Na_2CO_3	Na_2HPO_4	$\text{Na}_2\text{S}_2\text{O}_3$	$\text{Na}_2\text{S}_2\text{O}_5$
Na_2SO_3	Na_2SO_4	NaBr	NaCl	NaF	NaH_2PO_4
NaHCO_3	NaI	NaNO_2	NH_4Br	NH_4Cl	NH_4HCO_3
ZnCl_2	ZnSO_4	$(\text{BiO})_2\text{CO}_3$	$\text{Al}(\text{OH})_3$	AlPO_4	BaSO_4
$\text{Ca}(\text{OH})_2$	$\text{Ca}_3(\text{PO}_4)_2$	CaCO_3	CaSO_4	$\text{Mg}(\text{OH})_2$	MgCO_3
TiO_2	ZnO				

Code for the unknown analysis


Water soluble

Water insoluble

Code for the unknown analysis




Clicking the keys, you can type the six-character code of the analysis.



Settings menu




Here you can change some settings of the simulator.



Interaction with the objects¹

- The objects with which you can interact are highlighted on hover changing their color (orange or red) or size (like the book and the TV).







- You can interact with the objects only by clicking with the mouse or touching the screen (tap).
- The single click (tap) selects the object or starts an action.
- In the case of the trash beaker, double click (double tap) empties automatically all glassware containing something. Several video tutorials show this useful action.

Interaction with the objects²



- To use an object, it must be selected by a mouse click (tap). The selection is highlighted by lifting or moving of the object in the foreground (see reagents and analyte).



Interaction with the objects³

The objects can interact each other through actions given by:

- Selecting the object 1 and clicking the object 2 (e.g. selection of the stirrer and click the test tube to stir the content).
- Selecting the object, selecting the liquid reagent and clicking the reagent dropper (e.g. Selection of the test tube, selection of HCl 2M and click the reagent dropper).
- Selecting the object, selecting the solid reagent or the analyte and clicking the spatula (e.g. selection of the test tube, selection of NaNO₂ and click the spatula).



WARNING:

When an action is not permitted, the speech synthesis warn you and, if it is disabled, the action is just not completed without notice.

The glassware¹



The **glassware is always dirty** and needs to be rinsed with deionized water before the use. The simulator prevents the use of dirty glassware, avoiding the student from making mistakes in the analysis.

Generic cleaning procedure:

- Click the **glassware to clean** (e.g. a test tube) which will be then selected.
- Click the **wash bottle**. Some water will be transferred to the glassware.
- Click the **trash beaker**. The glassware content will be emptied.

Dropper cleaning procedure:

- Click an **empty test tube**. If you are not sure on its content, empty and wash it. The tube will be selected.
- Click the **wash bottle**. About 1 mL of water will be transferred to the tube.
- Click the **dropper** to select it.
- Click the **test tube with the water**. The dropper will suck up the tube content.



The glassware²



- Click the [trash beaker](#). The content will be released.

Ni-Cr wire cleaning procedure:

1. Turn on the [bunsen](#) setting the flame as oxidizing.
2. Put some [concentrated hydrochloric acid](#) in a test tube.
3. Click the [Ni-Cr wire](#) selecting it.
4. Click the [test tube with HCl](#). The wire will be soaked in HCl.
5. Click the [bunsen](#). The wire will be placed into the flame.
6. Repeat the steps from 3 to 5 until the flame does not change more the color (usually this operation must be repeated 3 times).



Transfer of substances

The simulator approximates the liquid and solid quantities according to the action and the size of the glassware on which you act:

Glassware	Liquid	Solid
Test tube	1 mL	One spatula tip
Beaker	5 mL	One spatula tip
Watch glass	1 mL	One spatula
Curcible	5 mL	One spatula
Paper	One drop	-

Therefore, if you want to add 2 mL of water to a test tube, you need to repeat the action twice because 1 mL is transferred at time.

There are some reactions in which the quantities play an important role for their occurrence as:

- **Borate test** (you must add methanol at least twice so that is more than sulfuric acid and can ignite).
- **Nitrate test** (you must add FeSO_4 at least twice so that the solution is enough concentrated to permit the formation of the brown ring).

Chemical reactions



LabSim is based on a database of about 3700 reactions that can occur between:

- Analyte + one reagent
- Analyte + two reagents
- Two reagents
- Reaction product + one reagent
- Reagent product + two reagents



The physical aspects influencing the result of a reaction are:

- **Heat** given by:
 - **water bath** (e.g. phosphate test, tests with TAA, dichromate paper, etc.);
 - **electric heater** (e.g. calcination with con KHSO_4).
- **Time.**



WARNING:

In order to occur correctly a reaction, the solution must be homogeneous after the addition of each reagent stirring the mixture.



is a web application developed in 2020-2021
by Alessandro Pedretti
All rights reserved

Alessandro Pedretti
Dipartimento di Scienze Farmaceutiche
Facoltà di Scienze del Farmaco
Università degli Studi di Milano
Via Luigi Mangiagalli, 25
I-20133 Milano - Italy
Tel. +39 02 503 19332
E-Mail: info@vegazz.net
WWW: <https://www.ddl.unimi.it/labsim>